Automatic Cleaning, Disinfection and Lubrication System

NSK

iCare+

iCare+

OPERATION MANUAL

Please read this Operation Manual carefully before use, and file for future reference

Powerful Partners®

"Congratulations for purchasing your new iCare+ device"

iCare+ is an innovative system whose function is to ensure the treatment of rotary dental devices (contra-angles, turbine) and dynamic dental devices. Its benefits are several:

- iCare+ insure cleaning-disinfection and lubrication of dynamic dental instruments
- iCare+ device allows the treatment of the inner and outer surfaces of different instruments (contra-angles, turbines...)
- Complete treatment in approx. 15 minutes
- Up to 4 instruments can be treated simultaneously and automatically
- The rotating transmission cleans, disinfects and lubricates even more effectively the contra-angles
- Validation and traceability insured by embedded software
- 3 types of treatment are available : full, partial, short

TABLE OF CONTENTS

CHAPTER 1: SAFETY INSTRUCTIONS	5
1.1. Intended use and user	5
1.1.1 Intended use	5
1.1.2 User	5
1.2. General safety measures	5
1.2.1. Adhesive safety labels	7
1.2.2. Position of safety labels	8
1.3. Individual safety measures	8
CHAPTER 2: DESCRIPTION OF iCare+ DEVICE	9
2.1. Contents of the sturdy box	9
2.2. Description of iCare+ device	10
CHAPTER 3: USE of iCare+ DEVICE	12
3.1. Setup	12
3.1.1. Installing the iCare+ device	12
3.1.2. Connection to the air system	12
3.1.3. Connecting to the electric network	14
3.2. Implementation of the products	15
3.2.1. Positioning of n.Clean and n.Cid bottles Implementation of the products	15
3.2.2. Filling the oil tank	17
3.3 Validation test	18
3.3.1. On site validation test	18
3.3.2. Periodical validation test	18
3.3.3. Routine test	18
3.4. Connection of the instruments	19
3.4.1. Treated instruments	19
3.4.2. Installation of the instruments	19
3.4.3. Removing the instruments	22
3.5. Programming of iCare+ device	24
3.5.1. The Control panel	24
3.5.2. Connecting the instruments	25
3.5.3. Selecting the program	25
3.5.4. Selecting the instruments	26
3.5.5. Selecting the type of instrument	26
3.5.6. Launching the cycle	28
CHAPTER 4: USER MAINTENANCE	30
4.1. General cleaning guidelines	30
4.2. General maintenance guidelines	30
4.2.1. Cleaning the treatment chamber	31

4.2.2. Cleaning the external parts _____	31
4.2.3. Emptying the used products of the drawer _____	31
4.2.4. Cleaning the drawer _____	32
4.2.5. Replacing the junction support blister _____	34
4.2.6. Replacing the O-ring joints of E-type connector _____	34
4.2.7. Replacing the fuse of the supply inlet _____	35
CHAPTER 5: ERROR MESSAGES FROM OPERATING LIGHT INDICATORS _____	36
CHAPTER 6: TRACEABILITY _____	37
CHAPTER 7: TROUBLE SHOOTING _____	37
CHAPTER 8: BREAK-DOWN GUIDE _____	38
8.1. Normal error _____	38
8.1.1. Liquid sensor _____	38
8.1.2. Oil sensor _____	38
8.1.3. Drawer sensor _____	38
8.1.4. Air sensor _____	39
8.1.5. Door sensor _____	40
8.2. Fatal error _____	40
CHAPTER 9: SPECIFICATIONS _____	41
9.1. Technical specifications: _____	41
9.2. Classification of the device _____	41
9.3. Symbol Mark _____	42
CHAPTER 10: ACCESSORIES AND PART LIST _____	43
CHAPTER 11: APPENDIXES _____	44
CE Certificate of Conformity _____	44
EMC INFORMATION _____	45
Reprocessing Instruction _____	47
WARRANTY _____	50

CHAPTER 1: SAFETY INSTRUCTIONS

1.1. Intended use and user

1.1.1 Intended use

The iCare+ is an automatic decontamination system intended to be used for disinfecting invasive medical devices in dentistry (soiled dynamic dental handpieces) prior to sterilization.

For this, iCare+ carries out four mandatory steps of a decontamination process, namely:

- Purge
- Cleaning
- Disinfection
- Lubrication

1.1.2 User

iCare+ device is a product intended to all qualified professionals from the dental field in medical typed environments such as private clinics, hospitals, dental offices or else university laboratories.

1.2. General safety measures

The user is advised to comply with the following instructions:

- The device must be used in accordance with the instructions given in this user's manual in terms of safety precaution and use of the system
- The power supply of the location where iCare+ device will be used has to comply with IEC requests, or with local regulations in force

The manufacturer, installer and importer are only responsible for the device's safety, reliability and performance if the aforementioned points have been respected.

Prior to a written request, the manufacturer is able to provide a technical description, diagrams and test instructions, which he considers to be useful in order for an authorized skilled technician to repair and maintain the repairable parts of the device.

IMPORTANT: Read this manual carefully before using your device. It contains important information for a correct use of your iCare+ device

<p>Caution symbol</p> 	<p>This symbol is designed to catch the reader's attention to any text containing vital safety information. You are asked to read closely and understand properly each of the sections concerned immediately before using your iCare+ device</p>
--	---

<p>Information symbol</p> 	<p>Text preceded by this symbol contains useful information for using your iCare+ device</p>
--	---

iCare+ device has been designed in accordance with current safety standards for similar devices. Despite this, it is essential that a number of important safety measures have to be applied in order to avoid any potential incident or even accident.

These main safety requests have to be taken into consideration by the user.

	<ul style="list-style-type: none"> • Use the iCare+ device only indoor • Do not install or use iCare+ device near a naked flame owing to the risk of explosion • Do not install or use iCare+ device in direct contact with sunlight • Ensure that the iCare + device is placed in a ventilated place • Store iCare+ device at temperature between 0°C and 50°C and humidity rate under 80% • Use iCare+ device in a temperate room (22-40°C) • Do not place n.Cid bottles near a heat source • Do not turn the device upside down • Place the device horizontally on a flat surface • Only use NSK products when operating iCare+ device (n.Clean as cleaning product, n.Cid as disinfection product, NSK Oil as lubrication product) • During the installation of iCare+ device, leave a space of 5 cm free on each side • Be sure that air supply pressure provided by iCare+ device is set between 5 and 6 bar • Do not remove the bottles (n.Clean or n.Cid) while iCare+ device is in use • If not using iCare+ device for a long period, switch of it and cut the air arrival • Only use NSK components for the maintenance of iCare+ device. The use of different components may damage the device
---	---

1.2.1. Adhesive safety labels

Adhesive safety labels have been placed for iCare+ device on important locations. These include storage instructions of the bottles and the device and information about the device's features.

	<p>The adhesive labels described below must be kept intact and should, if necessary, be replaced with the substitution labels provided with the accessories. In order to avoid damaging these labels, do not use abrasive products to clean your iCare+ device.</p>
---	---

 NSK iCare+ AUTOMATIC CLEANING - DISINFECTION LUBRICATION SYSTEM	REFERENCE: S102001 TYPE: I CARE + C2 SERIAL NUMBER: 0RYYYWNNNN    <div style="border: 1px solid black; padding: 2px; display: inline-block;">2015 </div>
SUPPLY VOLTAGE: 230VAC FREQUENCY: 50 Hz POWER INPUT: 50 VA FUSE RATING: 1.6A / 250V TIME LAG PRESSURE SUPPLY: 5 to 6 BAR	Manufactured by Canon Bretagne for: NSK Europe GmbH. Elly-Beinhorn Straße 8 D-65760 Eschborn, Germany  0086 Made in France CBO-E438-001

LA-4

LA-5

Figure 1: List of safety labels for iCare+ device

- LA-1: NSK OIL label (product for lubrication)
- LA-2: n.Clean label (product for cleaning)
- LA-3: n.Cid label (product for disinfection)
- LA-4: iCare+ device label with serial number
- LA-5: Warning label about risk that may be caused by the needles

1.2.2. Position of safety labels

The safety labels are glued on different parts of iCare+ device. The labels LA-1-2-3 are positioned on the internal surface of the top cover and as for LA-4 it is placed on the bottom panel of the unit housing.

The label LA-5 is positioned near the housing for the bottles of products.

1.3. Individual safety measures

All the safety measures referenced below must be followed closely in order to prevent exposure to risks of the user.

	<ul style="list-style-type: none"> • Never handle the power cord with wet hands; risk of electrocution • CAUTION: NSK does not warrant the use of a power cord other than that provided in the package. • Make sure the wall outlet has an earth connection with a differential (breaker) • Be careful not to inject water into the device; fire risk as a result of short circuits • Make sure to switch OFF the unit before any cleaning or maintenance action; risk of electrocution • Do not attempt to dismantle any part of the device. In case of breakdown, contact directly your retailer • In case of emission of smoke or a smell of burning, switch off immediately the iCare+ device and contact directly your retailer • Do not use the iCare+ with flammable gases
--	--

CHAPTER 2: DESCRIPTION OF iCare+ DEVICE

In this chapter, the main elements of iCare+ device and its features, useful for the user, will be described.

2.1. Contents of the sturdy box

Article n°	Description	Quantity
1	CD of installation with User Manual	1
2	Power cord	1
3	A 4/6 mm diameter tube for the pneumatic system (length: 2 m) with "T" junction connector	1
4	iCare+ device	1
5	n.Clean, n.Cid, lubrication oil (Products not included into iCare+ package. Must be ordered separately)	-

- At the delivery of the unit, please verify that the packaging has no damage during shipment.
- For shipment, always use the packaging of iCare+.
- Keep the packaging away from kids.

2.2. Description of iCare+ device

- Top cover for access to the housing of implementation of products: A
- Treatment chamber with connection of instruments: B
- Door for access to the treatment chamber: C
- Removable drawer for collecting the soiled products: D
- Control panel which enables the user to monitor the device: E

Figure 2: Front panel of iCare+ device

Figure 3: Rear panel of iCare+ device

ENGLISH

Figure 4: Top rear panel of iCare+ device + USB connection

- Power supply inlet: F
- Air filter: G
- Housing for the different products: H
- USB connection: I


Figure 5: Control Panel of the iCare+ device

- A: Cycle selection button A': Operating light indicators (LEDs) for cycle selected
- B: Ext w-o (or Spray) button for selection of the model of the handpiece(s)
- C: Instruments selection and Oil volume button (Short or Long mode)
- D: Start-up button. This button can also be used to stop the unit during cycle (with a long press).
- E: Cycle progress operating light indicators (LEDs)
- E: Warning operating light indicators (LEDs) for consumables state and for drawer, air, door state
- 1, 2, 3, 4: Indicator of the type of instrument (if the indicator is lit: instrument with external spray or without spray)

CHAPTER 3: USE of iCare+ DEVICE

CAUTION: The use of controls, settings or the implementation of procedures other than those specified below may expose users and patients to danger

3.1. Setup

3.1.1. Installing the iCare+ device

You should install the iCare+ device as requested below:

- The device should be positioned on a flat, solid and level surface. Do not level itself but the surface it is placed on. The maximum weight of the device with the maximum weight, fully loaded, in conditions of use is 16 kg
- The device should not be placed near a sink or any source that may splashed it
- The device should be placed in a room properly ventilated
- The device should be kept away from any source of heat

3.1.2. Connection to the air system

Insert the air supply tube in the air filter: J, fixed on the Rear panel's device. Be sure that the air tube is attached properly to the filter (see Figure 6). Connect the "T-shaped fast-on" connector to the surgery compressor's air supply and after that connect your dental unit.


Figure 6: Connection of the tube in the pneumatic system of iCare+ device

	<ul style="list-style-type: none"> • Following the norm EN ISO 7494-2, the air quality provided to iCare+ has to be dry, clean and free from bacteria and contamination. • The pressure of the air supply should be between 5 and 6 bar • The minimum air flow level provided to the unit is 50 L./mn. • It is necessary to follow the above mentioned parameters in order to avoid some dysfunctions or damage of the unit.
--	--


Figure 7: Connection of the iCare+ device to the air system

	<ul style="list-style-type: none"> • The pressure of the air supply must not be under 5 bar or above 6 bar. A pressure sensor detects if the air pressure is too low or too high and stops the machine. The operator is informed by an operating light indicator positioned on the Control Panel. • Do not fold or deform the air supply tube • If you experience some difficulties during the connection step, please contact your retailer • The air filter has an effective filter size of 5µm
--	---

3.1.3. Connecting to the electric network

Connect the device to the electric network by means of a power cord provided into the iCare+ package. This cable is a Class I cable.

Just insert the power cord in the socket: K, on the iCare+ device's panel.


Figure 8: Connection of iCare+ device to the electric network

CAUTION: NSK do not warrant the use of other power cord than the one provided with iCare+.

3.2. Implementation of the products

3.2.1. Positioning of n.Clean and n.Cid bottles Implementation of the products

Start opening the Top cover and take the bottles of n.Clean and n.Cid solutions supplied (see Figure 9).

Figure 9: n.Clean & n.Cid ready to be positioned on iCare+ device

	<p>CAUTION:</p> <ul style="list-style-type: none"> • Keep the products in a cool dry place and well ventilated. Not exposed to sunshine. • Avoid any contact with skin and eyes. • It is recommended to use protective glasses and gloves during handling of bottles or during emptying of waste drawer. • In case of contact with the eyes, rinse the eyes with water with the eyelids open for a sufficient length of time, and then seek ophthalmological advice immediately. • In case of contact with skin, wash off immediately with soap and water. In case of skin reactions, seek medical advice. • In case of inhalation, move affected person into fresh air. In case of irritation of the airways, seek medical advice. • If swallowed, rinse out mouth and then drink plenty of water. Do not induce vomiting. Seek medical advice.
--	--

ENGLISH

	<p>CAUTION: Do not forget to remove the metallic cap on each bottle, before starting the implementation of the consumables; you may damage the iCare+ device.</p>
--	--

Figure 10: Implementation of n.Clean bottle (500 mL)

Figure 11: Implementation of n.Cid bottle (500 mL)

Push the bottle of n.Clean solution into the correct housing (see Figure 10) and after that do the same thing with the n.Cid solution bottle (see Figure 11). To insert the bottles, a short pressure on them is enough. It is not necessary to push them strongly because they can be damaged.

	<p>For each bottle of solution or for the Oil, a coloured sticker is used to avoid confusions between the products during the implementation of them:</p> <ul style="list-style-type: none"> • For the n.Clean solution, a yellow sticker on its bottle and a yellow safety adhesive label on its housing (of the Top rear panel) are glued • For the n.Cid solution, a grey sticker on its bottle and a grey safety adhesive label on its housing (of the Top rear panel) are glued • The oil is directly poured into the tank (on the Top rear panel) on which a blue safety adhesive label is glued
--	---

	<p>Always use n.Clean and n.Cid bottles supplied by NSK Europe GmbH. The efficiency of the iCare+ device is optimized and certified for n.Clean and n.Cid products. Cleaning and disinfecting protocols have been calculated using these products. The use of any other product may lead to a failure or damage of the iCare+ device or of the instruments connected on it, and the efficiency of the device will not be ensured. Therefore, in case of non respect, the user will be considered as responsible.</p>
--	--

	<p>CAUTION: When the bottles are not inserted, the needles, if reachable, can cause a risk of injury</p>
--	---

3.2.2. Filling the oil tank


Start shaking well the maintenance Oil. Unscrew the tank's cap (see Figure 12) and pour the Oil in (see Figure 13). Please fill the tank up to 20 mm from the top edge of the tank. After filling, screw its cap.

Figure 12: Unscrewing the Oil tank's cap of iCare+ device

Figure 13: Filling the Oil tank of iCare+ device

- Always use maintenance NSK Oil for handpiece in the iCare+ device to ensure optimal operating conditions
- After filling the tank, be sure to have properly closed its cap

Once the tank is filled and the cap properly screwed, close the Top cover. After that, you can start connecting the different instruments.

3.3 Validation test

3.3.1. On site validation test

This test must be performed prior to first start-up. It consists of a series of tests in order to validate that the functioning parameters are well respected after the installation of iCare+.

3.3.2. Periodical validation test

We recommend to perform this periodical validation test every year, in order to validate that the performance of iCare+ is still optimum.

3.3.3. Routine test

A protein residual value determination using real contaminated instruments is part of the validation (chapter 3.3.1). NSK recommends a periodic, quarterly, protein residual value determination outside the revalidation cycles. The procedure is as follows:

- Attaching instruments from the practice circle after previous preparation in iCare+
(See the instructions for use chapter 12)
- Start the processing cycle
- Stop the processing cycle after the cleaning step has been completed and remove the instruments

Prior to the protein determination to be carried out according to the manufacturer, ensure that all fluids present in the instrument and on the instrument are completely removed (using compressed air and using a suitable adapter) in order to achieve a correct result. It is important to ensure that no cross-contamination occurs (e.g. by replacing the protection gloves for each instrument, etc.)

Afterwards, the instrument has to be completely reprocessed again in iCare +.

3.4. Connection of the instruments

3.4.1. Treated instruments

iCare+ is a device whose function is to ensure the cleaning and disinfection of dental rotary and dynamic instruments.

iCare+ device allows the cleaning and disinfection of internal and external surfaces of different instruments. The main instruments intended to be treated are usual handpieces: contra-angles and turbines.

The different types of instruments that can be treated by iCare+ device are listed in the table 1 (please refer to the [chapter 3.5.5](#). Selecting the type of instrument).

	<p>Caution: The instruments must be reprocessed following NSK and local authority requirements.</p> <ul style="list-style-type: none"> • External pre-cleaning and disinfection manually • Automated cleaning and disinfection process internal and external • Sterilization <p>For the manual of reprocessing refer to chapter 11</p>
---	--

ENGLISH

3.4.2. Installation of the instruments

With iCare+ device, it is possible to connect 2 kinds of instruments by means of 2 different connectors:

- Air turbine connector for turbines for dental.
- E-type connector for contra-angles for dental.

To connect a turbine with Midwest 4 holes(in accordance with the ISO 9168 standard) , start opening the door and connect the turbine directly to turbine connector of the iCare+ (see Figure 14) and screw the ring on tightly (see Figure 15). Finish this operation by closing the door. Take care of the position of the turbine connector.

	<p>Do not use the door as a support: do not put anything on the door when open.</p>
---	---

Figure 14: Connection of one turbine

Figure 15: Screwing the ring on the Midwest coupling

To connect turbine with other connection type, it is necessary to install first, the corresponding adaptor (Please refer to the list of available adaptors into the chapter 10).

Figure 16: Change of turbine adaptor

Insert the adaptor and screw the ring on firmly; check the positioning of the different nozzles.

ENGLISH

It is necessary to use a specific iCare+ air turbine adaptor. These are especially designed to deliver the right volume of product required for each instrument. NSK does not guarantee a good working of the device with any other coupling.

To connect a contra-angle or a hand piece on e-type connectors, in accordance with the ISO 3964 standard, start opening the door and connect the contra-angle on the E-type connector (see Figure 17). Finish this operation by closing the door.

Figure 17: Connection of a contra-angle

Figure 18: iCare+ version 1 with 1 air turbine and 3 E-type connectors

Place properly the contra-angle in the connector until you hear a “click”. In this way, the contra-angle will be well maintained during the different phases of the treatment.

CAUTION:

- Please make sure to remove any bur from the clutch
- Make sure also to remove any support/holder for snap-on caps
- If not the dirt and liquid inside the instrument could not come out

3.4.3. Removing the instruments

For removing an air turbine, in accordance with the ISO 9168 standard, follow the description below:

Once the cycle has finished, you will hear a "BIP" sound and the green operating light indicator ("END" LED) will be lit, start opening the door. After that, remove the turbine from the connection.

Removing a NSK turbine from its adaptor:

- Move up the retaining ring of the adaptor.
- Remove the turbine while maintaining the ring in upright position.

Retaining ring

Removing a Kavo turbine from its adaptor:

- Remove the turbine by drawing down the instrument.

The operation is identical to remove a turbine Sirona type of its corresponding adapter.

Removing a W&H turbine from its adaptor:

- Move up the retaining ring of the adaptor.
- Remove the turbine while maintaining the ring in upright position.

Removing a Bien-Air turbine from its adaptor:

- Press the button on the adapter.
- Remove the turbine while keeping the button pressed.

For removing a contra-angle, in accordance with the ISO 3964 standard, follow the description below:

Start opening the door. After that, remove the instrument in the axial direction of the E-type connector by pressing on the release button.

3.5. Programming of iCare+ device

Start putting the ON/OFF switch, localized on the Rear panel, on position I.

3.5.1. The Control panel

The Control panel (see Figure 19) is composed of a series of control buttons and of light operating indicators (LEDs).

The configuration of the device, the progress in the cleaning process and the status of the sensors are directly displayed on the front panel using 25 LEDs. Combinations of the LEDs provide additional information such as the progress in the copy of the internal memory into the USB key.

A total of 7 switches are available. One switch: (labelled "Cycle") is used to select one of the 3 types of cycle to run: Full cycle (Cleaning + Disinfection + Lubrication), partial cycle (Cleaning + Lubrication) or short cycle (lubrication only).

The "ext-WO" button enables the user to select the type of handpieces: a handpiece without spray or external spray, or by default with an inner spray.

4 switches enable to select/deselect the instruments to be processed and the volume of lubricant to be used.

The "Start" switch is used to start the programmed cycle or as an emergency stop.


Figure 19: Control Panel of the iCare+ device's Front panel

- A: Cycle selection button À: Operating light indicators (LEDs) for cycle selected
- B: Ext w-o (or Spray) button for selection of the model of the handpiece(s)
- C: Instruments selection and Oil volume button (Short or Long mode)
- D: Start-up button. This button can also be used to stop the unit during cycle (with a long press).
- E: Cycle progress operating light indicators (LEDs)
- F: Warning operating light indicators (LEDs) for consumables state and for drawer, air, door state
- 1, 2, 3, 4: Indicator of the type of instrument (if the indicator is lit: instrument with external spray or without spray)

3.5.2. Connecting the instruments

The instruments have to be connected in accordance to the description done in [chapter 3.4](#).

3.5.3. Selecting the program

The treatment of the instruments is achieved by a cycle composed of 3 different steps:

- 1 cleaning phase : 7 min (for 4 instruments)
- 1 disinfection phase : 7 min (for 4 instruments)
- 1 lubrication phase : 1.5 min (for 4 instruments)

The “CYCLE” button enables the operator to choose between the 3 maintenance cycles:

- Full maintenance cycle: Cleaning/Disinfection/Lubrication
- Partial maintenance cycle: Cleaning/Lubrication
- Short maintenance cycle: Lubrication only (similar function as for Care3plus)

Full mode

Partial mode

Short mode

Figure 20: Selection of the type of program for iCare+ device

A 4th pressure on the Cycle button manages the user to reset the choice of the cycle.

3.5.4. Selecting the instruments

By default, and as a measure of security, the iCare+ device dispenses the product to all the connectors. Here the operator has to select the instrument(s) and choose the suitable Oil volume that he wants to be delivered:

- Small drop corresponds to a normal volume (Short cycle)
- Large drop corresponds to a specific volume for certain contra-angles (Long cycle for surgical contra-angle or 1:5 contra-angle for example)
- 2 unlit LEDs indicate that the instrument is deselected and so that no product will be sent through it

	<p>The volumes and the sequences of injection of n.Clean and n.Cid products are directly programmed and monitored by the device itself, according to an accurately defined and validated protocol by a medical institute.</p>
---	---

By default, the Short cycle of oil is selected; press the button corresponding to the position where the instrument is connected, once to select the Long cycle (large drop LED lit) and twice to deselect the instrument.

3.5.5. Selecting the type of instrument

In parallel, the user will be able to choose the type of instrument.

The iCare+ device is designed to enable the treatment of 2 models of dental rotary instruments:

- Dental devices with only one drive channel
- Dental devices with 2 channels: W/A + drive channels

All the different instruments that can be treated with iCare+ device are mentioned in the following table:

	INSTRUMENT 1	INSTRUMENT 2	INSTRUMENT 3	INSTRUMENT 4
ONLY DRIVE CHANNEL	/	Contra-angle or handpiece with external spray or without spray		
W/A + DRIVE CHANNELS (2 CHANNELS)	Air turbine Air Scaler	Contra-angle or handpiece with internal spray		

Table 1: List of different treatable instruments per model

Standard setting of iCare+ is designed for internal spray instruments. To select for external or without spray instrument, it is necessary to push simultaneously the “ext/WO button” and the button corresponding to the position of the instrument without or external spray.

Button for selection of External/Without spray instrument

For example to change the model of instrument 2, the operator will have to press simultaneously the button "ext W/O" and the button 2; see on Figure 21 below. The LED "ext W/O", near button 2, is now lit to inform the operator that the instrument 2 is an instrument without spray or with external spray.


Figure 21: iCare+ device with lit W/O LEDs

In parallel of the function of external spray instrument, the unit can control if there is a mistake for the programming of the button SPRAY or if a channel is blocked.

For example, if the user selects an instrument as an external spray instrument (but in reality, the instrument has an internal spray), the unit is able to detect it.

As a reverse case, if the user connects an instrument with external spray without using the button "ext W/O" (to select external spray: "ext-w/o" LED), the unit is able to detect that the selection is wrong because this instrument has external spray.

To inform the user of these problems, the unit will stop and the Air LED will blink, with a "BIP" signal heard. In parallel, the concerning instrument is indicated by white LED blinking (oil indicator and/or ext-WO LEDs).

To restart the unit, first open the door, remove the instrument (to reduce the remaining pressure) reconnect the instrument, and then close the door.

If the problem happens again, remove the handpieces and push the button “Start”.

Check the instruments (in case of possible blocked channels).

3.5.6. Launching the cycle

Once the program is chosen, press START button to launch the cycle (see Figure 22).

	<ul style="list-style-type: none"> START button allows starting the programmed cycle for the selected instruments. If a problem happens, or that you wish to stop the unit during a cycle, START button allows also stopping the unit by a long press. The cycle is then stopped, the settings and selections are kept and the machine stays in standby. By pressing again START button, the cycle is relaunched from the beginning with the kept settings. If the door is open while the unit is in standby, the settings and selections are initialized.
--	--

ENGLISH

Figure 22: Launching the treatment cycle

	<p>If the user needs to stop the cycle during process, a long press on “START” button will stop the unit.</p>
--	---

A series of LEDs, positioned on the Control panel and on the treatment chamber, indicates the progress of the cycle:

- Yellow LEDs are used to provide information of cleaning phase
- White LEDs are used to provide information of disinfection phase
- Blue LEDs are used to provide information of lubrication phase

Once the cycle completed, green LEDs ("END" LED + treatment chamber LED) are lit and inform the user that the instruments are ready to be used again.

	<ul style="list-style-type: none"> • When the "Mode" light indicator (small/big amount of Oil) is switch off, the connection (1, 2, 3 or 4) is not selected • Just after purchasing iCare+ device, or after a long period of no use, run the device empty (without instruments)
--	---


Figure 23: Treatment cycle of iCare+ device

	<ul style="list-style-type: none"> • In order to facilitate the external cleaning, we recommend using disinfectant wipes (without protein fixation)for handling the instruments while connecting them to the unit and rub slightly them before connecting in case of important soils. It will help to liquefy the soils and stains located on the outside surfaces. • At the end of the cycle, open the door and remove the instrument with a disinfectant wipe (without protein fixation).
--	---

CHAPTER 4: USER MAINTENANCE

There are two different types of maintenance that have to be distinguished; one more global maintenance, that is regularly realized by the user and the specific maintenance performed by a skilled technician from NSK support.

Only skilled technician from NSK or people trained by NSK, are allowed to carry out specific maintenance.

4.1. General cleaning guidelines

	<ul style="list-style-type: none"> • CAUTION : All the cleaning operations should be performed after having disconnected the iCare+ device from the main power supply • Do not smoke in the room where the device is kept
---	---

	<p>You should clean the iCare+ device as requested below:</p> <ul style="list-style-type: none"> • Always use a soft cloth to clean the metallic and plastic parts of the device • Only use detergents that contain a very low level of alcohol (with no protein fixation agent) • Never attempt to clean areas, that are difficult to access, with sharp objects • Be careful when cleaning the Control panel and never use aggressive detergents
--	--

4.2. General maintenance guidelines

FREQUENCY / NUMBER OF CYCLES		OPERATION	PART NUMBER	DESCRIPTION
Every Week	50	Cleaning the treatment chamber	--	4.2.1.
Every Week	50	Cleaning the external parts	--	4.2.2.
2-3 times per week	15	Emptying the used products of the drawer	--	4.2.3.
Every Week	50	Cleaning the drawer	--	4.2.4.
quarterly		Replacing O-Rings an instrument holder	On demand	4.2.6.
1 year		Replacing the junction support blister	S103205	4.2.5.
When blown		Replacing the fuse(s) of the supply inlet	On demand	4.2.7.

Table 2: Frequency of maintenance for iCare+ device

4.2.1. Cleaning the treatment chamber

The treatment chamber can be easily cleaned, since the door is totally removable. Thanks to the 2 spring-loaded tabs, the door can be easily removable ensuring a sufficient access to the chamber for the cleaning step.

After that, if necessary disconnect and remove the respective instruments that are still implemented. Use a soft cloth or a sponge to clean the chamber and after that to rinse it. Proceed in the same way with the door.

We advise to use alcohol wipes for this part of the device. It will ensure to remove with efficacy the possible drops of stains.


Figure 24: Access to the treatment chamber of iCare+ device

4.2.2. Cleaning the external parts

Clean the external parts with a soft cloth, as mentioned in the part 4.1. Do not use scouring agents or highly abrasive products.

4.2.3. Emptying the used products of the drawer

The drawer should regularly be emptied. The emptying of the drawer is done after each use, at the end of the day.

In parallel when the drawer is full, the corresponding LED on the Control panel is lit.

Tray + Drawer

Figure 25: Removal of the drawer of iCare+ device

First start removing the drawer from its housing by sliding it backwards (see Figure 25) and secondly to dispose of the used products, undo the clips to remove the tray from the drawer (see Figure 26).


Figure 26: Removal of the tray clips and of the drawer's tray


Figure 27: Tray and drawer after removal

ENGLISH

	<p>All the products used for the treatment cycle of iCare+ device : Oil and especially the n.Clean and n.Cid solutions are formulated to be disposed of in the sewage disposal</p>
--	--

4.2.4. Cleaning the drawer

To avoid the creation of some residues or else to remove them, it is advised to clean regularly the drawer. According to the table 2 of part 4.1, the frequency for tray's cleaning should be weekly.

To realize this step you have to proceed in 2 steps: first, remove the clips on both sides of the tray and then remove the cover of the drawer (see Figure 27). After that, you can clean it by means of alcohol wipes.

4.2.5. Replacing the junction support blister

The blister allows the collecting of the possible drops of products during its change; it had to be replaced periodically.

As mentioned in the part 4.1. General cleaning guidelines, the frequency of replacement of the support junction blister is 1 year.

First of all open the top cover of iCare+ device to remove the soiled blister (see Figure 28), and then put directly the new one on the support junction (by being sure that this last is properly sealed by means of the new blister).

	<p>CAUTION: When removing the soiled blister, do not empty out the used products on the iCare+ device, on the surface of its and on the floor.</p>
---	---


Figure 28: Replacement of the junction support blister

4.2.6. Replacing the O-ring joints of E-type connector

Press the O-ring joint and remove it from the E-type connector (see Figure 29).

Thread the new O-ring joints in the corresponding grooves.

Figure 29: O-ring joints of E-type connector

4.2.7. Replacing the fuse of the supply inlet

As mentioned in the part 4.1, the fuse(s) has (have) to be changed when blown.

Press down simultaneously and pull the two tabs at the end of the fuse housing in order to unlock it. Remove and replace the fuses according to their rating. Replace the fuses in the fuse box and insert them in their housing. Push to lock (see Figure 30).

Figure 30: Fuse box of iCare+ device

Remark:

The fuses should be replaced functions of the country of use; the specifications for each type of fuse are summarized on the following table.


NATURE OF FUSE	VOLTAGE	AMPERE RATING
EUROPE	230 VAC	T1.6 AH 250V

For all the other maintenance operations of your iCare+ device, please contact immediately your NSK retailer.

CHAPTER 5: ERROR MESSAGES FROM OPERATING LIGHT INDICATORS

Description of operating light indicators:

The different indicators and the Control panel are useful to inform the operator in real time about the operating conditions of the iCare+ device.

- A. “Cle” indicator informs the operator that n.Clean product level is too low. iCare+ device stops. It is necessary to replace the bottle by a new one in order to restart the device. Remove the used bottle with a firm upward movement. To install a new bottle refer to paragraph 3.2.1.1.
- B. “Cid” indicator informs the operator that n.Cid product level is too low. iCare+ device stops. It is necessary to replace the bottle by a new one in order to restart the device. Remove the used bottle with a firm upward movement. To install a new bottle, refer to paragraph 3.2.1.1
- C. “Oil” indicator informs the operator that the lubrication Oil product level is too low. iCare+ device stops. It is necessary to fill the oil tank. To fill the oil tank, refer to paragraph 3.2.1.2.

After changing a bottle (n.Clean or n.Cid) or filling the Oil tank, press “Start” to prime the circuit.

- D. “Drawer” indicator informs the operator about 2 things:
 - Either the drawer is not properly positioned (or even absent). In this case, check that the drawer is present and correctly inserted.
 - Or the drawer is full and might overflow. In this case, remove the drawer and empty it.
 In both cases, the iCare+ is unable to start the cycle.
- E. This light indicator informs the operator that there is a problem with the air pressure (too low air pressure value, no air and too high air pressure value). Air pressure plays an important role in the cycle efficacy; iCare+ device cannot start the cycle until the problem is not solved.
- F. This light indicator informs the operator that the door is open or not properly closed. Close the door to start any cycle.

iCare+ is equipped with a microprocessor that checks in real time every internal component. If any problem is detected on one of the component, iCare+ device stops immediately

CHAPTER 6: TRACEABILITY

The cleaning and disinfection phases have to be performed in accordance with the different current recommendations. It is therefore necessary to enable the validation of the handpieces treatment done by the iCare+ device. That's why traceability software must be used. The software "SoftCare+" was created to ensure the traceability monitoring of the cycles realized by the device.

Please refer to specific user manual of annex software destined to visualize cycles information.

CHAPTER 7: TROUBLE SHOOTING

PROBLEM(S)	POSSIBLE CAUSE(S)	SOLUTION(S)
The device does not switch on	Is the supply cord properly connected to mains?	Connect the supply cord to mains
	Is the supply cord properly connected to the iCare+ device?	Insert properly the cord plug in the device's main connector
	Is the switch positioned on the "ON" position?	Turn the main switch to "ON" position
	Are the fuses working?	Change the fuses with the same model
The device does not clean and/or lubricate	Are the mode light indicators lit?	Push on one of the selection cycles buttons
	Is the door opened? (door LED lit)	Close the door properly
	Are the n.Clean or n.Cid level LEDs lit?	Change the corresponding bottle and prime
	Is the Oil level LED lit?	Fill the Oil tank
The device does not clean and/or lubricate correctly	Is the air supply tube properly connected to the device?	Connect properly the air supply tube to the device's connector
	Is the air supply pressure set between 5 and 6 bar?	Set the air supply pressure between 5 and 6 bar
	Is Short mode selected?	Select Long mode
	Is the instrument properly connected to the connector?	Plug the instrument properly

CHAPTER 8: BREAK-DOWN GUIDE

The other problems, not previously listed in the chapter 5 (Error messages from operating light indicators) and in the chapter 7 (Possible failures), are explained in the following chapter.

8.1. Normal error

8.1.1. Liquid sensor

Corresponding LED is lit and a “BIP” signal can be heard. Please change the bottle and push “START” button to fill the product line.

If the bottle is not empty, the sensor could have detected a bubble into the product line. Please push “START” button to fill the product line.

8.1.2. Oil sensor

Corresponding LED is lit and a “BIP” signal can be heard. Please open the oil tank and fill it with maintenance oil.

8.1.3. Drawer sensor

Corresponding LED is lit and a “BIP” signal can be heard. Remove the drawer and empty it. Or else, make sure that the drawer is in appropriate position.

8.1.4. Air sensor

With this air sensor, there are 2 LED signals. Corresponding LED is lit constantly with “BIP” signal or corresponding LED is blinking with “BIP” signal.

- If corresponding LED is constantly lit and a “BIP” signal can be heard.

The air pressure is lower than 5 bar or higher than 6 bar.

To restart the unit, first, check the air pressure at the inlet of the machine. It is necessary to have between 5-6 bar at the inlet. Once this is done, press “START” button to reset the error.

- If “Air” LED is blinking with the LEDs around the buttons 1-2-3-4 (oil quantity selection + ext-WO) and a “BIP” signal can be heard.

This measurement of dynamic air pressure is done to detect if channels are blocked.

This means 3 possibilities of alert:

- Tube blocked inside the unit
- Channel blocked inside the instrument
- Wrong programming in case of use of ext-WO button for external spray instrument

Tube blocked inside the unit:


The unit can detect if one of the tube inside the machine is blocked.

Blocked channel/External spray function:

The unit can detect if one of the channel of the instrument is blocked (spray channels or drive channels)

The iCare+ indicates the concerned instrument by blinking the oil quantity and ext-WO LEDs next to the corresponding button (1-2-3-4)

In parallel of the function of external spray instrument, the unit can control if there is a mistake for the programming of the ext-WO button.

Button for selection of External spray /without spray instrument

For example, if the user select an instrument as an external spray (but in reality, the instrument has an internal spray), the unit is able to detect it.

As a reverse case, if the user connects an instrument with external spray without using the button SPRAY (to select external spray: “ext-w/o” LED), the unit is able to detect that the selection is wrong because this instrument has external spray.

Whatever the reason of the alert is (blocked channel or wrong selection), to restart the unit, it is first necessary to check the reason of the problem. Then, open the door, disconnect the instruments and reconnect them. Finish by closing the door.

Once this is done, push “START” button to reset the error. The “Air” LED should be switch off.

Check the instrument, if one of the channels is not blocked, and program the right selection (internal or external spray instrument) by using the ext-WO button + the number of the instrument.

ENGLISH

8.1.5. Door sensor

Corresponding red LED “Door” is lit (No “BIP” signal).

Please, make sure that the door is properly closed.

8.2. Fatal error

The 6 red LEDs are blinking and 3 “BIP” sounds are heard.

If this case happens, please contact directly your retailer.

CHAPTER 9: SPECIFICATIONS

9.1. Technical specifications:

ELECTRICAL SUPPLY	Power supply voltage: 230 VAC Input frequency: 50 Hz
iCare+ DEVICE Use environment (working temperature): Transportation and storage environment (storage temperature): Use environment (humidity rate): Transportation and storage environment (humidity rate): Maximal altitude of use: Use environment (atmospheric pressure range): Transportation and storage environment (atmospheric pressure range): Mains voltage: Power input : Pollution level Dimensions overall: Maximum space required: Weight empty: Maximum weight, fully loaded, in conditions of use: Maximum heat output: Maximum noise level: Air pressure level : Air filtration :	At temperate room 22-40°C 0-50°C Specific maximal humidity of 80% for temperatures to 31°C, with a linear decrease to a specific humidity of 50% at 40°C 10-80% 2000 m 700-1060 hPa 500-1060 hPa Variation shall not be over +/- 10% of the nominal voltage 50 VA II H: 405 mm / W: 355 mm / D: 400 mm H: 455 mm / W: 455 / D: 450 mm 14 kg 16 kg Not significant < 60 dBA >5 bar and <6 bar 5µm
PRODUCTS	500 mL packaged in a PEHD bottle to be pierced 500 mL packaged in a PEHD bottle to be pierced Tank of 500 mL to be filled
OTHERS	Fully micro-processor driven and controlled (embedded software)/Programming of the iCare+ device thanks to an interface (Control panel)

9.2. Classification of the device

The iCare+ device belongs to the II b class in accordance with the rule 15 of the annex IX mentioned in the council directive 93/42/EEC regarding medical devices.

9.3. Symbol Mark

	Mandatory conformance mark on many products placed on the market in the European Economic Area
	Conforms to CE European Directive of “Medical equipment directive 93/42/EEC.” (PATENT PENDING)
	Manufacturer.
	Follow the waste of electric and electronic equipment (WEEE) Directive (2002/96/EC) for product and accessory disposal
	Consult operation instructions.
	Caution, Refer to attached instructions
	For indoor use only.
	Use by
	Keep Dry.
	Lot number
	Indication of the range of humidity level
	Indication of pressure range
	Indication of temperatures range
	Indication of the year of manufacturing
	Protective gloves recommended
	Protective glasses recommended

CHAPTER 10: ACCESSORIES AND PART LIST

The accessories and the spare parts for iCare+ device are respectively listed in the tables 3 and 4.

PRODUCT DESCRIPTION	REFERENCE
n.Clean, 500 mL PEHD bottle (6pcs) (not included into iCare+ package. Must be ordered separately)	ACL600
n.Cid, 500 mL PEHD bottle (6pcs) (not included into iCare+ package. Must be ordered separately)	ACD600
Maintenance Oil, 1 L PEHD bottle (not included into iCare+ package. Must be ordered separately)	Z016117
iCare+ PTL Adaptor (not included into iCare+ package. Must be ordered separately)	Z1127010
iCare+ KV Adaptor (not included into iCare+ package. Must be ordered separately)	Z1127011
iCare+ SR Adaptor (not included into iCare+ package. Must be ordered separately)	Z1127012
iCare+ WH Adaptor (not included into iCare+ package. Must be ordered separately)	Z1127013
iCare+ BA Adaptor (not included into iCare+ package. Must be ordered separately)	Z1127014
4/6 mm diameter tube (1 piece included into iCare+ package)	On demand
Power Cord (1 piece included into iCare+ package)	On demand
CD of operation manual (1 piece included into iCare+ package)	On demand

Table 3: List of accessories for iCare+ device

PRODUCT DESCRIPTION	REFERENCE
E-type O-ring joint (Set of 4 pieces)	Y900-185
Blister pack	S103205

Table 4: List of spare parts for iCare+ device

PTL Adaptor refers to adaptor(s) to use for the connection of NSK[®] air turbines for PTL-CL-4HV-T during the treatment cycle.

KV Adaptor refers to adaptor(s) to use for the connection of KAVO[®] handpieces during the treatment cycle.

SR Adaptor refers to adaptor(s) to use for the connection of SIRONA[®] handpieces during the treatment cycle.

WH Adaptor refers to adaptor(s) to use for the connection of W&H[®] handpieces during the treatment cycle.

BA Adaptor refers to adaptor(s) to use for the connection of BIEN-AIR[®] handpieces during the treatment cycle.

CHAPTER 11: APPENDIXES

CE Certificate of Conformity

The iCare+ device complies with the following standards:

- 61010-1 and 61010-2-040 standards: Safety requirements for electrical equipment for measurement, control and laboratory use:
 - Part 1: General requirements,
 - Part 2-040: Particular requirements for sterilizers and washer-disinfectors used to treat medical materials

This product has been granted CE Marking in accordance with European directive (EN/93/42/CE) that is applicable to laboratory and more precisely to medical devices.

The number 0086 placed under the CE label refers to the notified body which has validated the product's conformity and so issued the certificate of conformity.

The head office responsible for the product is the following one:

NSK Europe GmbH
Elly-Beinhorn-Straße-8
D-65760 Eschborn
GERMANY

EMC INFORMATION

The iCare+ device complies with the following standards:

- IEC 61326-1 standard: Electrical equipment for measurement, control and laboratory use-EMC requirements; Part 1: General requirements

Guidance and manufacturer's declaration - electromagnetics emissions		
iCare+ device is intended for use in the electromagnetic environment specified below. The user of iCare+ should assure that is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - Guidance
RF Emissions EN 55011	Group 1	iCare+ should not cause any interference in nearby electronic equipment.
RF Emissions EN 55011	Class B	iCare+ is suitable for use in all establishments, including domestic establishments and those directly connected to public low-voltage power supply network that supply network that supplies buildings used for domestic purposes.
Harmonics Emissions EN 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions EN 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetics immunity			
iCare+ device is intended for use in the electromagnetic environment specified below. The user of iCare+ should assure that is used in such an environment.			
Immunity test	EN 61326-1 Test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharges (ESD) EN 61000-4-2	+/- 4kV contact +/- 4kV air	+/- 4kV contact +/- 4kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst EN 61000-4-4	+/- 1kV for power supply lines	+/- 1kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge EN 61000-4-5	+/- 0.5kV line(s) to line(s) +/- 1kV line(s) to earth	+/- 0.5kV line(s) to line(s) +/- 1kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines EN 61000-4-11	0% Ut during 10ms 0% Ut during 20ms 70% Ut during 200ms 0% Ut during 5 sec.	0% Ut during 10ms 0% Ut during 20ms 70% Ut during 200ms 0% Ut during 5 sec.	Mains power quality should be that of a typical commercial or hospital environment. If the user of iCare+ requires continued operation during power mains interruptions, it is recommended that the iCare+ be powered from an uninterruptible power supply.
Magnetic field EN 61000-4-8	2 V/m	2 V/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Guidance and manufacturer's declaration - electromagnetic immunity			
iCare+ device is intended for use in the electromagnetic environment specified below. The user of iCare+ should assure that is used in such an environment.			
Immunity test	EN 61326-1 Test level	Compliance level	Electromagnetic environment - Guidance
Conducted RF EN 61000-4-6	3V eff 150kHz to 80MHz	3V rms	Portable and mobile RF communications equipment should be used no closer to any part of iCare+ including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ de 80MHz à 800MHz $d = 2.3\sqrt{P}$ de 800MHz à 2.7GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter
Radiated RF EN 61000-4-3	3V /m 80MHz to 2.7GHz	3V/m	

			<p>manufacturer and is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters are determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b). Interference may occur in the vicinity of equipment marked with the following symbol :</p> 
--	--	--	---

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobiles radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which iCare+ is used exceeds the applicable RF compliance level above, the iCare+ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the iCare+.

b Over the frequency range 150kHz to 80MHz, field strengths should be less than 3V/m.

Cables and accessories	Maximum length	Complies with
AC cord	2m	RF emissions Harmonic emissions 2 EN 61000-3-2 Voltage fluctuations / flicker emission 3 EN 61000-3-3 Electrostatics discharges (ESD) 2 EN 61000-4-2 Electric fast transient / burst 4 EN 61000-4-4 Surge 5 EN 61000-4-5 Voltage dips, short interruptions and voltage variations on power supply input lines 11 EN 61000-4-11 Magnetic fields 8 EN 61000-4-8 Conducted RF 4-6 EN 61000-4-6 Radiated RF 3 EN 61000-4-3

ENGLISH

Recommended separation distance between portables and mobile RF communications equipment and the iCare+

iCare+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The user of iCare+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment and the iCare+ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150kHz to 80MHz d=1.2VP	80MHz to 800MHz d=1.2VP	800MHz to 2.5GHz d=2.3VP
0.001	0.12	0.12	0.23
0.01	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitter rated at a maximum output power not listed above, the recommended separation distance « d » in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where « P » is the maximum output power rating of the transmitter in watts (W) according to the transmitters manufacturer.

NOTE 1 At 80 Mhz and 800Mhz, the separation distance for the higher frequency is applicable.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Reprocessing Instruction

Recommendation by RKI

Medical devices that are contaminated with pathogens can be a source of infection for humans.

Thus, the use of such medical devices requires previous reprocessing, for which defined requirements have been defined by Robert Koch Institute (RKI).

Before reprocessing, it is necessary:

- to verify the product compatibility with the reprocessing method to be used
- to check how often the reprocessing has to be performed
- to perform a biological risk assessment in order to determine the classification of the medical device

Once performed, reprocessing shall comprise the following individual steps to be performed according to the classification of the instruments:


- Proper preparation
- Inspection
- Cleaning
- Disinfection
- Maintenance and repair
- Functional testing
- Identification
- Packaging
- Sterilization

This document will only address Preparation, Cleaning and Disinfection with regards to intended use "iCare +".

1. Preparation of staff / protective measures against possible infections




- For your safety: Always wear protective gloves (1), mouth and nose protection as well as glasses to minimize the risk of infection during machine conditioning
- Consider national regulations and guidelines (e.g. ISO 17664, RKI, etc.)

Note: Steps 2 to 5 must be performed within 10 minutes after the end of treatment. If this is not possible, the instruments must be pre-treated with the non-protein-fixing cleaning solution IC100 (from Alpro Medical) according to the manufacturer's instructions.

2. Preparation of instruments

- Allow air and / or water to run on the instrument after treatment of a patient for at least 20 seconds (2)
- Immediately wipe off the instruments afterwards with an approved disinfectant (for example, MinutenWipes / Alpro) (do not spray on!) (3)
- Hand-/Contra-angled Instruments:
 1. Remove the instrument from the micromotor (4)
 2. Remove the drill from the chuck system (e.g. with a pair of tweezers) and immerse it in a disinfecting solution (5)
- Turbines:
 1. Disconnect the turbine from the coupling (or disconnect fixed connection of the turbine from the hose) (4)
 2. Remove the drill from the clamping system (e.g. with a pair of tweezers) and immerse it in a disinfecting solution (5)

Never immerse NSK handpieces & turbines in disinfecting solutions or in ultrasonic cleaning devices!


3. Contamination-proof transport

to the place of preparation

4. Pre-cleaning outside

- Pre-clean hand & angle pieces and turbines under running cold water using a soft brush until all visible soiling is removed (6)

5. Cleaning (internal & external), disinfection (internal & external) and subsequent oil treatment with NSK iCare +


- Lock the instruments to be processed on the adapters (7) of the iCare + device provided for this purpose and close the door
- Carry out the necessary settings on the device as per instructions for use and start the processing cycle
- Do not interrupt the cycle and remove the attached instruments only when the iCare + device signals the end of the processing cycle.

Only use n.clean (cleaning solution) and n.cid (disinfection solution) (8) as well as care oil from NSK to fill the iCare + device

- Remove the cleaned, disinfected and well-groomed instruments from the device after the completion of the processing cycle and dry them with a paper towel (9). If necessary, drain any excess oil or blow or drain the instruments for a few seconds before applying them to the patient

9. Visual Check


- Check the instrument for cleanliness and intactness

10. Sterilization

Packed sterilization (DIN EN ISO 13060) according to the autoclave manufacturers' load pattern following these parameters:

Autoclave B-Class:

At least 3 min. At 134 ° C alternatively: At least 20 min. At 121 ° C

WARRANTY

The manufacturer will warrant the quality of the product for one year after you purchased it, against manufacturing errors and defects in equipment under conditions of normal installation and use and of after-sales monitoring. NSK reserves the right to analyse and determine the cause of any problem. The parts considered as consumables are not covered by the manufacturer's warranty.

Note, however, that if you did not observe what is written in this operation manual or for any consumables or recommended liquids, the warranty will not apply.

**NSK Europe GmbH**

Elly-Beinhorn-Strasse 8,
65760 Eschborn,
Germany
Tel: +49 (0) 6196 77606-0
Fax: +49 (0) 6196 77606-29
www.nsk-europe.de/
info@nsk-europe.de

NSK France

32, rue de Lisbonne
75008 Paris,
France
Tel: +33 (0)1 56 79 59 80
Fax: +33 (0)1 56 79 59 81
<http://www.nsk.fr>
info@nsk.fr

NSK Dental Spain SA

C/ Módena, 43
El Soho - Európolis
28232 Las Rozas,
Madrid,
Spain
Tel: +34-91-626-61-28
Fax: +34-91-626-61-32
info@nsk-spain.es
<http://www.nsk-spain.es>

NSK United Kingdom Ltd

Office 5, Gateway 1000,
Arlington Business Park,
Whittle Way, Stevenage,
Hertfordshire, SG1 2 FP
Tel: 0800 634 1909
Tel Intl: +44 (0)1438 310 670
Fax: 0800 634 1910
Fax Intl: +44 (0)1438 310 671
<http://www.nsk-uk.com>

Specifications are subject to change without notice

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© Springer-Verlag 2012